DYNASUN

USER MANUAL
BEDIENUNGSANLEITUNG

F307

TABLE OF CONTENTS

- 1** GETTING STARTED
INTRODUCTION
ERSTE SCHRITTE
INIZIO
- 2** COPYRIGHT NOTICE
COPYRIGHT
HINWEIS ZUM URHEBERRECHT
COPYRIGHT
- 3** Safety instructions
Sicherheitshinweise
Consignes de sécurité
Istruzioni di sicurezza
- 4** OPERATION
FONCTIONNEMENT
BETRIEB
MONTAGGIO
- 5** TECHNICAL SUPPORT
SUPPORT TECHNIQUE
TECHNISCHE UNTERSTÜTZUNG
- 6** WARRANTY INFORMATION

1.
GETTING STARTED
INTRODUCTION
ERSTE SCHRITTE
INIZIO

1.1 English
Welcome

Thank you for choosing DynaSun! We are delighted that you've selected our product, designed with state-of-the-art technology, premium materials, and expert craftsmanship.

The innovative design of your DynaSun ensures professional-grade performance, making it a reliable companion for your demanding work on set. As a high-quality tool, it empowers creative professionals like you to achieve outstanding results effortlessly..

We hope you thoroughly enjoy using your DynaSun. To ensure maximum safety and optimal performance, please take the time to carefully read the operating manual.

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COPYRIGHT

2.1 English

This document contains proprietary information belonging to Confidence Europe – DynaSun Division. The content is protected by copyright and may not be reproduced, in whole or in part, without prior written consent from Confidence Europe.

While every effort has been made to ensure the information provided is accurate and up-to-date as of the publication date, no guarantees are made or implied regarding its completeness, accuracy, or freedom from errors. Confidence Europe reserves the right to make changes at any time for the purpose of product improvement.

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3.2 English

Safety Instructions for the Ergonomic Four-Leg Walking Stick

1. **Pre-Use Check:**
 - **Visual Inspection:** Check the walking stick for cracks, dents, or other damage before each use. A damaged walking stick should not be used.
 - **Feet and Rubber Tips:** Ensure the feet and rubber tips are in good condition to avoid slipping.
2. **Proper Use:**
 - **Height Adjustment:** Adjust the height of the walking stick to match your body height, ensuring the grip is at hip level when standing.
 - **Usage:** Hold the walking stick with one hand while using the other hand for additional support.
3. **Walking with the Stick:**
 - **Movement:** Move the stick in front of you and alternate placing your feet on the ground. Use the stick for support by pressing it down on its stable legs.
 - **Steps:** Take small, controlled steps to maintain balance.
4. **Safety in Daily Life:**
 - **On Slippery Surfaces:** Be cautious on slippery or uneven surfaces. Avoid walking on smooth floors without appropriate shoes.
 - **Stairs:** Use a handrail when going up or down stairs and only use the stick for additional support.

Maintenance of the Ergonomic Four-Leg Walking Stick

1. **Cleaning:**

- **Regular Cleaning:** Clean the walking stick regularly with a damp cloth and mild detergent. Avoid harsh chemicals that could damage the material.
- **Feet and Rubber Tips:** Regularly check the feet and rubber tips for wear and clean them to ensure good grip.

2. **Inspection:**

- **Regular Checks:** Regularly inspect the walking stick for loose parts or signs of wear. Replace damaged parts promptly.
- **Grip:** Ensure the grip is secure and stable. Replace it if there are signs of wear.

3. **Storage:**

- **Safe Storage:** Store the walking stick in a dry, secure place to avoid damage. Avoid placing it near heat sources or direct sunlight.

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English

Using the Ergonomic Four-Leg Walking Stick

1. **Preparation:**
 - **Height Adjustment:** Ensure the walking stick is set to the correct height before use.
 - **Getting Ready:** Stand upright before using the stick. Hold it with your dominant hand.
2. **Walking with the Stick:**
 - **Maintain Balance:** Keep the stick close to your body and place it on the ground while lifting your other foot to move forward.
 - **Step by Step:** Move the stick forward, followed by one foot, and repeat at a steady pace.
3. **Sitting and Standing:**
 - **Sitting:** Use the stick for stabilization as you sit down. Hold it with one hand while using the other to support yourself on the chair.
 - **Standing Up:** Gently push yourself up with the stick and the chair while leaning your body forward.

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Support technique

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5.1 English

If you require technical support for your product, please first contact the retailer where you purchased the device. In most cases, the retailer can provide quick and competent assistance with your questions or issues.

If the retailer is unable to assist you or recommends contacting our support team directly, our technical support team will be happy to help. Please use the contact information provided below to reach us.

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6.1 English

We want to hear from you!

Confidence Europe – DynaSun Division provides a twelve (12) month warranty from the date of purchase by the end consumer. This warranty covers defects in materials and workmanship in the mechanics, control electronics, and power supply under normal use and proper maintenance. Proof of purchase is required to make a warranty claim. If no proof of purchase is provided, the shipping date to the first purchaser (retail or distributor) will be used to determine the warranty period.

The warranty will be void if the product is not handled with care, especially in the case of damage caused by electrostatic discharge, adverse temperature and humidity conditions, or misuse. Confidence Europe – DynaSun Division reserves the right, at its discretion, to repair or replace products that are found to be defective within the warranty period and that are returned to Confidence Europe – DynaSun Division, F.O.B. the place of origin.

This warranty is exclusive, and there are no other written or oral, express or implied, warranties. In particular, Confidence Europe – DynaSun Division disclaims any implied warranty of merchantability or fitness for a particular purpose to the extent permitted by law.

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Confidence Europe GmbH